

C-flex AP-Vantage™

Radiolucent Head Positioner

The C-Flex AP-Vantage is the most versatile and easiest to use radiolucent head positioner on the market. It connects to any OR table, and produces the cleanest artifact-free cervical images possible.

Ideal for use in all cervical procedures including:

- Posterior Cervical
- Anterior Cervical
- Odontoid Screws

- A-70110 C-Flex AP-Vantage (with C-Flex System)
- A-70111 C-Flex AP-Vantage (without C-Flex System)

Specifications

- Fits Accessory Rail Types: US, EU, UK, Denyer, Japanese (specify rail type)
- Fits OR Tables with Outside Rail Width 20" - 25 1/2" (508 mm - 648 mm)
- Device Weight = 20 lbs. (9 kg)
- 500 lbs. (227 kg) Patient Weight Limit
- Works with the Mayfield Skull Clamp and Mayfield Horseshoe
- Includes:
 - OR Table Adaptor
 - C-Flex Traction Bar
 - C-Flex 45° Head Module
 - C-Flex AP-Vantage Adjustable Leg

C-flex Range of Motion Specifications

- Lateral Adjustment Range 6"
- Height Adjustment Range 28"
- Extension Adjustment Range 9"
- Lateral Tilt + 15°

C-flex Head Positioning System Components

- C-Flex Base Unit
- Starburst Adaptor
- Prone Mask Head Module
- Flat Plate Module
- Comfort Mask Disposables
- C-Flex Clean Cape Disposables

Allen Medical Systems, A Hill-Rom Company
One Post Office Square, Acton, MA 01720
(800) 433-5774

www.allenmedical.com

UNITED STATES

Toll-free (800) 433-5774
Phone (978) 263-7727
Fax (978) 263-8846

INTERNATIONAL (OUTSIDE US, CANADA)

Phone +44 (0) 870 112 0769
Fax +44 (0) 1530 561 174

MAYFIELD is a registered trademark of SM USA, Inc.
Jackson Tables are not affiliated with or manufactured by Allen® Medical Systems.

© 2008 Allen Medical Systems, A Hill-Rom Company. All Rights Reserved.
D-770396-A1 October 30, 2008

ALL PRODUCTS ARE CLASS 1 MEDICAL DEVICES

allen



C-flex AP-Vantage™

Radiolucent Head Positioner

allen


Artifact-Free Imaging...



Odontoid Screw




Concord Position



Jackson Table Imaging Capabilities On Any OR Table

Just Squeeze & Position...



C-flex Head Positioner: Cervical Flexion / Extension



Anterior



Posterior



Surgeon's View



- Optional Traction Bar
- C-flex Head Positioner
- OR Table Adaptor
- Optional Adjustable Leg



Fits Any OR Table

C-flex AP-Vantage Overview



C-flex Accessories

- 
Prone Mask
- 
Comfort Mask Disposable
- 
Flat Plate
- 
Foam Block Disposable
- 
Mayfield Clamp