

FlashPak®

Mobile Sterilization Container System

The Symmetry Riley FlashPak® is a hospital-proven closed container system engineered for safe transport of flashed instruments in a mobile sterile field. The FlashPak is appropriate for use with gravity and pre-vacuum cycles.

The FlashPak meets current guidelines for flash sterilization. Effective sterilization is achieved in five or 10-minute flash cycles. Sterilized instruments can be stored up to 24 hours inside the sealed system.

The completely sealed system prevents contamination during transport giving dependable protection you require.

The system accepts standard wire baskets and various trays to accommodate many instrument types. The transparent lid allows quick identification of instruments. The advanced plastic polymer construction provides guaranteed durability.

Contact your local authorized distributor for more information.



Protection, in a *Flash*

FLASHPAK® MOBILE STERILIZATION SYSTEM

Product Descriptions



- **FlashPak System**

System comes complete with Base, Wire Basket and Lid.

Small FlashPak Part #: 9020
11.5 x 11.2 x 7.3"
(29.2 x 28.5 x 18.5cm)

Medium FlashPak Part #: 9030
11.2 x 20.7 x 8.4"
(28.5 x 52.9 x 21.3cm)

Large FlashPak Part #: 9040
11 x 24.7 x 9"
(27.9 x 62.7 x 22.9cm)

Extra Large FlashPak Part #: 9050
13 x 24.7 x 9"
(33 x 62.7 x 22.9cm)



- **FlashPak Replacement Parts**


Wire Baskets

Small Wire Basket for 9020: 9020-08
Medium Wire Basket for 9030: 9030-08
Large Wire Basket for 9040: 9040-08
Extra Large Wire Basket for 9050: 9050-08



Additional Replacement Parts

Refurbish Kit: 9020-09 (2), 9020-11 (2): 9000-09
Valve for 9020,9030,9040,9050: 9020-09
Pressure Balance Element for all Bases: 9020-11
Balance Element Gasket for all Bases: 9020-12



- **FlashPak Accessories**

FlashPak Delivery Cart (all sizes): 9111
Autoclave Gloves (one pair): 9000-14

Silicone Mats

Silicone Pin Mat for 9020-08: 9020-04
Silicone Pin Mat for 9030-08: 9030-04
Silicone Pin Mat for 9040-08: 9040-04
Silicone Pin Mat for 9050-08: 9050-04

 **SYMMETRY Medical^{inc}**

27 Wrights Landing, Auburn, ME 04210, USA
T: +1.207.786.2775 Fax: +1.207.786.5553
www.symmetrymedical.com

Cat. No. 8002-00 Rev. A 2M707
Printed in USA. ©2007 Symmetry Medical Inc. All rights reserved.

 **SURGICO**
MEDICAL & SURGICAL LTD
P.O. Box 300 180 Albany
Freephone: 0800 400 333
Freefax: 0800 900 222