


# MICROPAK<sup>®</sup>

MICROSURGICAL INSTRUMENT TRAY



## Instrument protection in a complete range of sizes.

Riley MicroPak<sup>®</sup> instrument trays are designed for use with a wide variety of microsurgical and other delicate instruments. These durable, lightweight plastic trays are engineered for maximized protection and ease of handling, with transparent lids for quick instrument identification. Made of rugged, state-of-the-art plastic polymers, the Riley MicroPak is compatible for use with the majority of sterilization methods currently available.

**Riley**  
MEDICAL

[www.RileyMed.com](http://www.RileyMed.com)



# MICROPAK<sup>®</sup>

## MICROSURGICAL INSTRUMENT TRAY

### Maximized Protection for Microsurgical Instruments

- **Protects** - helps minimize costly repairs of your most delicate surgical instruments
- **Secure** - instruments are held safely in place with non-slip mats or silicone holding bars
- **Flexible** - MicroPak is available in a full range of sizes to accommodate various delicate instruments
- **Durable** - engineered with highly resilient plastic polymers to withstand the rigors of repeated use and sterilization
- **Convenient** - transparent lids allows easy identification of instruments during transport, sterilization or storage

Cat#	Description	Dimensions*		
		Width inches/(mm)	Length inches/(mm)	Depth inches/(mm)
0100	Base & Lid	1.6 (40.6)	2.6 (66.0)	0.5 (12.7)
0101	Base, Lid & Mat	1.6 (40.6)	2.6 (66.0)	0.5 (12.7)
0111	Base, Lid & Mat	2.0 (50.8)	6.0 (152.4)	0.8 (19.1)
0112	Base, Lid & Mat	2.5 (63.5)	9.7 (246.4)	1.0 (25.4)
0126	Base, Lid & Silicone Bars(1 diamond knife)	2.0 (50.8)	6.0 (152.4)	0.8 (19.1)
0127	Base, Lid & Silicone Bars(2 diamond knife)	2.0 (50.8)	6.0 (152.4)	0.8 (19.1)
0128	Base, Lid & Silicone Bars(3 diamond knife)	3.0 (76.2)	6.0 (152.4)	0.8 (19.1)
0211	Base, Lid & Silicone Mat	3.0 (76.2)	6.0 (152.4)	0.8 (19.1)
1010	Base & Lid	6.0 (152.4)	10.0 (254)	0.8 (19.1)
1011	Base, Lid & Silicone Mat	6.0 (152.4)	10.0 (254)	0.8 (19.1)
1015	Base, Lid & Mat	5.5 (139.7)	9.5 (241.3)	2.5 (63.5)
1020	Base & Lid	6.0 (152.4)	10.0 (254)	1.5 (38.1)
1021	Base, Lid, Insert Tray & 2 Mats	6.0 (152.4)	10.0 (254)	1.5 (38.1)
1023	Base, Lid, Insert Tray & 1 Insert Mat	6.0 (152.4)	10.0 (254)	1.5 (38.1)
1025	Base, Lid & Mat	6.0 (152.4)	10.0 (254)	1.5 (38.1)
2011	Base, Lid & Mat	10.0 (254)	15.0 (381)	0.8 (19.1)
2020	Base & Lid	10.0 (254)	15.0 (381)	1.5 (38.1)
2021	Base, Lid, Insert Tray & 2 Mats	10.0 (254)	15.0 (381)	1.5 (38.1)
2023	Base, Lid, Insert Tray & 1 Insert Mat	10.0 (254)	15.0 (381)	1.5 (38.1)
2024	Base, Lid, Insert Tray & 2 Mats	10.0 (254)	15.0 (381)	1.5 (38.1)
2025	Base, Lid & Mat	10.0 (254)	15.0 (381)	1.5 (38.1)

- Notes:**
1. Dimensions are inside dimensions and are denoted as W x L x D in inches and millimeters
  2. Standard colors are Blue Base, Amber Lid, Amber Insert, Blue Mat

**Warranty:** Riley Medical, Inc. warrants the MicroPak<sup>®</sup> microsurgical trays against failure due to manufacturing defects. If the product(s) fail to function due to manufacturing defects within one year from the date of purchase, Riley Medical, Inc. will repair or replace the product(s) at no charge, to the customer upon prepaid return of the product(s) to Riley Medical, Inc. This warranty, which is given expressly and in lieu of all other warranties expressed or implied of merchantability and fitness for particular purpose, constitutes the only warranty made by Riley Medical, Inc.

Manufacturing a full line of engineered trays that protect and organize your instruments

**FLASHPAK<sup>®</sup>**  
Mobile Sterilization  
Container System

**MICROPAK<sup>®</sup>**  
Microsurgical Instrument Trays

**MULTIPAK<sup>®</sup>**  
Instrument Organizing Trays

**SCOPEPAK<sup>®</sup>**  
Endoscope Trays

**TRANSPAK<sup>®</sup>**  
Instrument Soaking  
& Transport Trays



**Riley Medical, Inc.**  
27 Wrights Landing, Auburn, Maine 04210-9915 USA  
1-800-245-3300 • (207) 786-2775 • FAX (207) 786-5553  
www.rileymed.com