

Shoulder Access System™

Arm Positioner

The Allen Shoulder Access System is the smallest, lightweight arm positioner on the market – for any shoulder procedure in the Beach Chair position. It attaches over the drape and provides effortless positioning of the arm perioperatively - at a fraction of the price of the Spider and the McConnell. Additionally, the disposable is more intuitive and competitively priced.

Ideal for use in all shoulder and arthroscopic procedures including:

- Rotator Cuff Repair
- Bankart Repair
- SLAP Lesion Repair
- Capsular Release

Specifications

- Fits Accessory Rail Types: US, EU, UK, DEN (specify rail type)
- Device weighs 6 lbs. (2.7 kg)
- 500 lb (227 kg) maximum patient weight
- A-92000 Shoulder Access System (US Rail)
- A-92000E Shoulder Access System (EU Rail)
- A-92000UK Shoulder Access System (UK Rail)
- A-92000D Shoulder Access System (Denyer Rail)
- A-92001 Shoulder Access System - Sterile Disposables (Case of 10)

Optional Accessories

The Allen Shoulder Access System works well with the Allen Lift-Assist Beach Chair or the Allen Table-Powered Beach Chair.



Allen Lift-Assist
Beach Chair (A-91000)



Allen Table-Powered
Beach Chair (A-90000)

Allen Medical Systems, A Hill-Rom Company
One Post Office Square, Acton, MA 01720
(800) 433-5774

www.allenmedical.com

UNITED STATES

Toll-free (800) 433-5774
Phone (978) 263-7727
Fax (978) 263-8846

INTERNATIONAL (OUTSIDE US, CANADA)

Phone +44 (0) 870 112 0769
Fax +44 (0) 1530 561 174

MAYFIELD is a registered trademark of SM USA, Inc.
Jackson Tables are not affiliated with or manufactured by Allen® Medical Systems.

© 2008 Allen Medical Systems, A Hill-Rom Company. All Rights Reserved.
D-770395-A1 October 29, 2008

ALL PRODUCTS ARE CLASS 1 MEDICAL DEVICES

Allen

Allen



Shoulder Access System™

Arm Positioner

Spend Your Time On The Procedure, Not The Positioning...



Unlimited surgical site access in the sterile field.



Lightweight / Portable



Effortless Positioning: Complete range of motion for all shoulder procedures.



Disposable: Intuitive and easy to use



Connects OVER the drape in the sterile field.

The Shoulder Access System Is Fully Autoclavable.