



SCOPEPAK[®]

ENDOSCOPE TRAYS

Protects endoscopes during transport, storage and sterilization

The Riley ScopePak[®] offers engineered protection for the handling and storage of costly endoscopes. Available in a variety of sizes to fit a broad range of endoscopes, each tray features silicones holding bars to secure the telescope in place, and a transparent lid for easy identification. Made of rugged, state-of-the-art plastic polymers, the Riley ScopePak is compatible for use with the majority of sterilization methods currently available.

Riley
MEDICAL

www.RileyMed.com



SCOPEPAK®

ENDOSCOPE TRAYS

**RILEY
MEDICAL**

*Manufacturing a full line
of engineered trays
that protect and organize
your instruments*

FLASHPAK®
Mobile Sterilization
Container System

MICROPAK®
Microsurgical Instrument Trays

MULTIPAK®
Instrument Organizing Trays

SCOPEPAK®
Endoscope Trays

TRANSPAK®
Instrument Soaking
& Transport Trays

Engineered Protection for Your Delicate Endoscopes

- **Proven Protection** - helps prevent damage to endoscopes, minimizing expensive repairs
- **Secure** - endoscopes are fixed in place with unique silicone holding bars
- **Flexible** - adjustable silicone holding bars accommodate a variety of endoscope sizes
- **Durable** - engineered with tough, time-tested, highly resilient plastic polymers to withstand the rigors of repeated use and sterilization
- **Convenient** - transparent lids allows easy identification of endoscopes

Catalog#	Description	Dimensions*		
		Width inches/(mm)	Length inches/(mm)	Depth inches/(mm)
9312	ScopePak, small w/Scope Holder Sets ¹	3.0 (76.2)	11.1 (281.9)	1.3 (33.0)
9317	ScopePak, medium w/Scope Holder Sets ²	3.0 (76.2)	16.1 (408.9)	1.3 (33.0)
9321	ScopePak, large w/Scope Holder Sets ²	3.0 (76.2)	20.1 (510.5)	1.3 (33.0)
9326	ScopePak, extra large w/Scope Holder Sets ³	4.0 (101.6)	26.0 (660.4)	2.8 (71.1)
9331	ScopePak, extra long w/Scope Holer Sets ³	4.7 (119.4)	30.3 (769.6)	2.8 (71.1)

Notes: 1. Includes 1 ea. #5350 and 1 ea. #5351 silicone scope holder
2. Includes 2 ea. #5350 and 1 ea. #5351 silicone scope holder
3. Includes 4 ea. #5353 and 2 ea. #5352 silicone scope holder
4. Dimensions are inside dimensions

Standard Colors: Base: blue; Lid: amber; Silicone Bars: blue

Warranty: Riley Medical, Inc. warrants the ScopePak® endoscope trays against failure due to manufacturing defects. If the product(s) fail to function due to manufacturing defect within **one year** from the date of purchase, Riley Medical, Inc. will repair or replace the product(s) at no charge, to the customer upon prepaid return of the product(s) to Riley Medical, Inc. This warranty, which is given expressly and in lieu of all other warranties expressed or implied of merchantability and fitness for particular purpose, constitutes the only warranty made by Riley Medical, Inc.



Riley Medical, Inc.
27 Wrights Landing, Auburn, Maine 04210-9915 USA
1-800-245-3300 • (207) 786-2775 • FAX (207) 786-5553
www.rileymed.com