

TRANSPAK®

INSTRUMENT SOAKING & TRANSPORT TRAY



Now there's a safer, more secure way to transport contaminated instruments.

The Riley TransPak® was engineered to protect personnel from the hazards of contaminated instruments after surgical use. This strong, lightweight plastic soaking system is designed for minimized spillage and splashing, helping protect hands and eyes from possible infection during transport to cleaning facilities. Its removable insert tray permits the safe transfer of contaminated instruments with no hand contact; the TransPak itself can be transferred directly to a washer/decontaminator.

Riley
MEDICAL

www.RileyMed.com



TRANSPAK®

INSTRUMENT SOAKING & TRANSPORT TRAY

Transport Contaminated Instruments with Confidence

Safe

- minimized spillage and splashing protects personnel
- removable insert trays allow transfer of contaminated instruments with no hands contact
- container clearly identified in Biohazard Red

Versatile

- made of special plastic polymers tough enough to be put directly in a washer/decontaminator
- used with cold sterilizing solutions
- engineered to withstand the rigors of repeat use

Convenient

- removable insert tray designed to accept Riley MultiPak® and most other instrument trays
- base and insert completely nest for ease of storage

*Manufacturing a full line
of engineered trays
that protect and organize
your instruments*

FLASHPAK®
*Mobile Sterilization
Container System*

MICROPAK®
Microsurgical Instrument Trays

MULTIPAK®
Instrument Organizing Trays

SCOPEPAK®
Endoscope Trays

TRANSPAK®
*Instrument Soaking
& Transport Trays*

Catalog#	Description	Dimensions*		
		Width inches/(mm)	Length inches/(mm)	Depth inches/(mm)
8600	TransPak w/ Base, Lid & Insert	15.0 (381.0)	25.8 (654.1)	6.5 (165.1)
8600-00	Base & Lid			
8600-01	Base only			
8600-02	Lid only			
8600-03	Insert Only			

Color: *Biohazard Red*

Notes: *Dimensions are inside dimensions.
Not recommended for steam sterilization.*



Warranty: Riley Medical, Inc. warrants the TransPak® against failure due to manufacturing defects. If the Product(s) fail to function due to manufacturing defect within 90 days from the date of purchase, Riley Medical, Inc. will repair or replace the product(s) at no charge, to the customer upon prepaid return of the product(s) to Riley Medical, Inc. This warranty, which is given expressly and in lieu of all other warranties expressed or implied of merchantability and fitness for particular purpose, constitutes the only warranty made by Riley Medical, Inc.



Riley Medical, Inc.
27 Wrights Landing, Auburn, Maine 04210-9915 USA
1-800-245-3300 • (207) 786-2775 • FAX (207) 786-5553
www.rileymed.com