

# When Every Second Counts, Expect More From Your Solidifier.



Contain Control  
**ViraSorb®**  
Super Solidifier  
Dispose

# ViraSorb® Super Solidifier... When every second counts, expect more from your solidifier!



## Faster

Gels in seconds ... up to 10 times faster than other solidifiers!

## More Effective

Evenly and completely gels canister contents without stirring or agitating. Does not cause fluid expansion.

## Stays Solid

Facilitating safer handling and transport to disposal site by guarding against spills or leaks.

## Concentrated Formula

Powerful formula surpasses the competition – 25 grams solidifies up to 1,000cc of fluids.

## Ideal Product Chemistry

Glutaraldehyde-free, chlorine-free and phenol-free formula.

## Safe

Non-toxic.  
Has no offensive odor.  
Safer to handle than fluids.  
Can be safely incinerated.

## Complies with Department of Transportation Packaging Requirements

ViraSorb® does not revert back to a fluid state, eliminating spillage, and leakage, qualifying it as *Solid Regulated Medical Waste*.

**Controls** fluid by quickly penetrating the canister contents to gel up to ten times faster than other solidifiers without stirring or agitation.

**Contains** fluid by staying solid, facilitating safer handling and transport of medical waste.

**Dispose** of in accordance with your facility's protocol, state and local regulations as *Solid Regulated Medical Waste*.

Ask about the WaterBoom™ Quiet Floor Suction Strip, another outstanding fluid management product from Colby!



For more information or to place an order, CALL 800-969-3718, FAX 800-897-2192

## Ordering Information

Item No.	Description	Case Quantity
1005	ViraSorb® Super Solidifier Solidifies up to 500cc of fluid	200 Dissolvable Pouches
1000	ViraSorb® Super Solidifier Solidifies up to 1,200cc of fluid	96 Single-Use Bottles
1001	ViraSorb® Super Solidifier Solidifies up to 1,500cc of fluid	96 Single-Use Bottles
1002	ViraSorb® Super Solidifier Solidifies up to 2,000cc of fluid	36 Single-Use Bottles
1003	ViraSorb® Super Solidifier Solidifies up to 3,000cc of fluid	36 Single-Use Bottles
1010	ViraSorb® Super Solidifier Solidifies up to 10,000cc of fluid	12 Multi-Use Bottles
1020	ViraSorb® Super Solidifier Solidifies up to 20,000cc of fluid	12 Multi-Use Bottles

Pour nozzles are included with Multi-Use Bottles. Pre-load pouches into suction canister for a closed system.



## Fluid Waste Management Products

**COLBY**  
manufacturing  
corporation

1016 Branagan Drive  
Tullytown, PA 19007  
[www.colbymfg.com](http://www.colbymfg.com)

ViraSorb® is a registered trademark of Colby Manufacturing Corporation. Colby Manufacturing Corporation is registered with the FDA and is certified to the ISO-9001:2000 Standard. Made in the USA. Printed in the USA. 1541 Rev E