

**You Worry About
What's On The Table;
And Let Us Worry
About What's Under It.**



**The WaterBug[®]
Quiet Floor Suction Device**

You worry about what's on the table; and let us worry about what's under it.



Wet floors in the OR are a hazard. The WaterBug® Quiet Floor Suction Device efficiently suctions fluid waste to reduce risks and maintain a dry operating room floor.

Low Profile
12' of flexible tubing exiting from the side reduces hazards.

Easy to Maneuver
Glides easily and will not stick to floors so you can easily reposition it.

Thorough Fluid Pick-up
Fluid is captured from the top via drain holes as well as from the bottom via suction and wicking action.

Complies with OSHA's Bloodborne Pathogen Standard.

Quiet Operation
Absorbent material ensures quiet wicking operation and easy mobility.

Control fluids generated during very wet procedures by easily maneuvering the WaterBug with a push of your foot to the areas affected by fluids. Its raised design allows for easy repositioning, controlling fluids over a large floor area.

Contain unlimited amounts of fluids from the floor with the WaterBug's thorough fluid pick-up design from the top via drain holes as well as from the bottom.

Dispose of the WaterBug and potentially hazardous waste in accordance with your facility's protocol for a quick and easy room clean-up.

The WaterBug® Quiet Floor Suction Device is a convenient, effective means of dealing with infectious waste fluid generated during very wet procedures. It suctions unlimited amounts of irrigants, blood and other body fluids from the floor. Using the WaterBug helps reduce the likelihood of injuries associated with wet floors. In addition, it helps protect healthcare professionals and housekeeping personnel from contact with contaminated substances and aids in your efforts to comply with OSHA's Bloodborne Pathogen Standard.

If solidifying, ViraSorb® Super Solidifier may be used to contain and solidify fluids to help your facility comply with the Department of Transportation's packaging requirements.



For more information or to place an order, CALL 800-969-3718, FAX 800-897-2192

Ordering Information

Item No.	Description	Case Quantity
90030	WaterBug® Quiet Floor Suction Device	30/case
90020	WaterBug® Quiet Floor Suction Device	20/case
90010	WaterBug® Quiet Floor Suction Device	10/case

Fluid Waste Management Products
COLBY
manufacturing
corporation
1016 Branagan Drive
Tullytown, PA 19007
www.colbymfg.com

WaterBug® is a registered trademark of Colby Manufacturing Corporation. Colby Manufacturing Corporation is registered with the FDA and is certified to the ISO 9001:2000 Standard. Made in the USA. Printed in the USA. 1103 Rev E