


An Innovative Stirrup For Today's Modern OR™

# The KYRA COMFORT™ Stirrup

## Introductory Offer \$9995 + Gst



**1x Pair OR Kyra Comfort™ or Velcro Stirrups & 1x Pair Intergrated Clamps. Additional Stirrup Carts are available for \$995 + Gst**

- **NO MORE LOST CLAMPS MAKES LIFE EASIER IN TODAY'S MODERN OR:**  
The Integrated KYRA Secure-Release™ Clamp solves the "lost clamp problem". Easy to mount, easy to release, the Integrated KYRA Secure-Release™ Clamp makes life easier and less hectic for clinical staff.
- **MORE PRESSURE REDUCING PADDING TO ACCOMMODATE LARGER LEGS OF TO DAY'S PATIENTS:**  
Kyra Medical's unique stirrup design allows 1/2" more padding for pressure relief than traditional stirrups. Unlike traditional boot pads, Kyra's pad protects the medial and lateral aspect of the ankle.
- **ELIMINATES DANGEROUS PINCH HAZARDS:**  
KYRA Comfort Stirrup's curved support rods eliminate dangerous pinch hazards that traditional stirrups present but without the use of hard to clean guards. The unique "blade design" also eliminates potential pinch hazards present when mounting traditional stirrups onto table rails.
- **BETTER INFECTION CONTROL:**  
Optional KYRA Secure-Lok™ Closure Straps eliminate Velcro® BRAND straps and are designed to improve infection control.